



screw tubes & tube accessories

▶ 15 & 50mL tubes

▶ pestles, openers & locks

15 & 50 mL centrifuge tubes

Sterile 15 and 50mL screw cap tubes fit all popular centrifuges. They can all be easily labelled with standard marker pens in the frosted or printed marking area. As only a ¼ turn of the cap is needed to safely and securely close and seal, one-handed operation is straightforward. All tubes are available in both loose, bagged configuration as well as in handy recyclable card racks.




15 and 50mL tubes



pestles, assorted



cap openers



cap locks

Item #	Description	Unit pack	Case pack
--------	-------------	-----------	-----------

Polypropylene tubes: sterile

▶ rated to 12,000 x G ▶ operating temperature range -90°C to +120°C

2935-SP-50400 15mL graduated

2935-SP-51100 15mL graduated

2935-SP-50100 50mL graduated

2935-SP-61400 50mL graduated

See Supplement Page S8

Polystyrene tubes (not autoclavable): sterile DISCONTINUED

▶ rated to 1,500 x G ▶ operating temperature range -40°C to +90°C

2837-SS-50400 15mL graduated bag of 50 10 units

2837-SS-51100 15mL graduated 2 racks of 25 10 units

Pestles

The DNase, RNase and DNA free pestles may be used to homogenize cells or resuspend precipitates such as protein or DNA pellets. The pestle & tube combination pack enables you to benefit from using the tube designed and manufactured to give a firm match and therefore achieve maximum efficiency with the pestle.

0.5mL pestle

1004-39 pestle only, assorted 10 per unit 40 units

1014-39 pestle with 0.5mL tube (SSI#1110-00) 10 of each 20 units

1.5mL pestle

1005-39 pestle only, assorted 10 per unit 40 units

1015-39 pestle with 1.5mL tube (SSI#1210-00) 10 of each 20 units

Openers and locks

Cap openers ease fatigue when opening and closing large numbers of microtubes. These useful, brightly colored tools also help in the safe opening of glass ampules.

Cap locks are designed to securely seal SSI tubes against the rigors of freezing and boiling. The large, easy-to-grip tab permits safe insertion and removal from boiling water and freezers. An assortment of bright colors together with a writing area allow easy identification or organization into types.

cap openers

1007-09 cap openers, assorted 100 10 units

cap locks

1002-39 for 0.5mL tubes, assorted 100 10 units

1003-39 for 1.5 or 2.0mL tubes, assorted 100 10 units